

## IRIDEX EndoProbe® Devices Verified with the Constellation\* Vision System

The following EndoProbe devices have been verified for use with the Constellation\* Vision System.



1212 Terra Bella Avenue  
Mountain View, CA  
94043-1824 USA

Phone: 650.940.4700  
Fax: 650.940.4710  
Web: www.igidex.com

EndoProbe® Models					
Type	Description	19.5 gauge	20 gauge	23 gauge	25 gauge
Standard	Straight		10562	14390	13920
	Angled 45°		10547		
Aspirating	Active Straight		14000		
	Active Angled		14010		
	Passive Fluted		11473		
Illuminating	Bayonet Straight		14420		
	Bayonet Angled 30°		14410		
	BriteLight™ Straight	13900		14540	14490
	BriteLight™ Angled 30°	14020			
	BriteLight™ Angled 45°	13930			
	Stepped			14545	14560
Adjustable & Intuitive	Finger Adjust		14572F	14573F	14574F
	Thumb Adjust		14572T	14573T	14574T
Stepped	Angled 45°		14030	14400	14120
RFID EndoProbe® Models					
Standard	Straight		65692	65716	65713
	Angled 45°		65695		
Aspirating	Active Straight		65752		
	Active Angled		65755		
Illuminating	Bayonet Straight		65728		
	Bayonet Angled 30°		65731		
	BriteLight™ Straight			65707	65704
	Stepped			65710	65950
	Angled 45°		65743	65698	65701

To use with an IRIDEX EndoProbe® device, please follow the steps below:

1. Plug an IRIDEX probe into laser port 1 or 2
2. Press **UNRECOGNIZE** on the monitor
3. Select **ENDO** for probe
4. If using an illuminating probe, select **ILLUMINATOR** Window
5. Select proper **GAUGE** and associated **WIDE FIELD OPTION**
6. Select **SURGERY** then **LASER** for main operating window

For additional information or to order IRIDEX EndoProbe® devices, please call 800-388-4747.

Sincerely,

Paul Hardiman  
VP of Regulatory